





Probe IGH/BCL6, DY Translocation Dual Color FISH Probe
Product No. TR054

IGH

Probe Size Approximately 1550kb
Fluorophore FITC
Location 14q32
Also known as IGD1; IGH@; IGHJ; IGHV; IGHD@; IGHJ@; IGHV@; IGH.1@; IGHDY1

BCL6

Probe Size Approximately 790kb
Fluorophore TexRed
Location 3q27.3
Also known as BCL5; LAZ3; BCL6A; ZNF51; ZBTB27



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.