




Probe API2/MALT1, DY Translocation Dual Color FISH Probe  
Product No. TR057

#### API2

Probe Size Approximately 1350kb  
Fluorophore FITC  
Location 11q21  
Also known as AIP1; BIRC3; MIHC; CIAP2; HAIP1; HIAP1; IAP-1; MALT2; RNF49; c-IAP2

#### MALT1

Probe Size Approximately 890kb  
Fluorophore TexRed  
Location 18q21  
Also known as MLT; MLT1; IMD12; PCASP1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.