





Probe ETV6/NTRK3, DY Translocation Dual Color FISH Probe
Product No. TR014

ETV6

Probe Size Approximately 910kb
Fluorophore TexRed
Location 12p13
Also known as TEL; THC5; TEL/ABL

NTRK3

Probe Size Approximately 1200kb
Fluorophore FITC
Location 15q25
Also known as TRKC; GP145-TrkC; gp145(trkC)



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.